# **High Frequency Lock-In Amplifier**

SR844 — 200 MHz, dual phase RF lock-in amplifier



- · 25 kHz to 200 MHz frequency range
- · 80 dB dynamic reserve
- Time constants from 100 µs to 30 ks (6, 12, 18 or 24 dB/oct rolloff)
- "No Time Constant" mode (10 μs to 20 μs update rate)
- · Auto gain, phase, reserve and offset
- Two 16-bit DACs and ADCs
- · Internal or external reference
- · GPIB and RS-232 interfaces

# SR844 200 MHz Lock-In Amplifier

The SR844 is the widest bandwidth lock-in amplifier available. Unlike simple down-converters, no additional instruments are required. And unlike analog lock-ins, no manual frequency range switching is necessary.

The SR844 provides uncompromised performance with a frequency range of 25 kHz to 200 MHz and up to 80 dB of drift-free dynamic reserve. It includes the many features, ease of operation and programmability that you've come to expect from SRS DSP lock-in amplifiers.

## **Digital Technology**

The SR844 uses the same advanced DSP technology found in the SR850, SR830 and SR810 lock-in amplifiers. DSP offers many advantages over analog instruments—high dynamic reserve, low zero-drift, accurate RF phase shifts and orthogonality, and digital output filtering.

## **Signal Input**

The SR844 has both 50  $\Omega$  and 1 M $\Omega$  inputs. The 1 M $\Omega$  input is used with high source impedances at low frequencies, or with a standard 10× scope probe. The 50  $\Omega$  input provides the best RF signal matching. Up to 60 dB of RF attenuation or 20 dB of RF gain can be selected in 20 dB increments. Full-scale sensitivities range from 1 Vrms (+13 dBm) to 100 nVrms (-127 dBm). Gain allocation can be optimized to provide up to 80 dB of dynamic reserve.



phone: (408)744-9040 www.thinkSRS.com

#### Reference

The SR844 offers both external and internal reference operation. In both cases, the entire 25 kHz to 200 MHz frequency range is covered without any manual range selection. The external reference input has an auto-threshold feature which locks to sine, square or pulsed signals. The internal reference is digitally synthesized and is adjustable with 3-digit frequency resolution.

Harmonic detection of the 2F component is available for both internal and external reference modes.

A reference output (1.0 Vpp square wave into 50  $\Omega$ ), which is phase synchronous with the lock-in reference, is available in both external and internal mode.

#### **Output Filters**

Time constants from 100  $\mu$ s to 30 ks can be selected with a choice of 6, 12, 18 or 24 dB/oct rolloff. For high bandwidth, real-time outputs, the filtering can be by-passed entirely. In this "No Filter" mode, the effective time constant is about 30  $\mu$ s with the analog outputs updating every 10 to 20  $\mu$ s.

#### **Ease of Operation**

The SR844 is easy to use. All instrument functions are set from the front-panel keypad, and the knob is used to quickly adjust parameters. Up to nine different instrument configurations can be stored in non-volatile memory for fast, reliable instrument setup. Standard RS-232 and GPIB (IEEE-488.2) interfaces provide connections to your data acquisition systems.

#### **Useful Features**

Auto-functions allow parameters that are frequently adjusted to be set automatically. Sensitivity, dynamic reserve, phase and offset are each quickly optimized with a simple key stroke.

The offset and expand features are useful for evaluating small fluctuations in your signal. The input is nulled with the auto-offset function and output expand increases the resolution by up to  $100\times$ .

Ratio mode is used to normalize the signal to an externally applied analog voltage. It is useful to eliminate the effect of source intensity fluctuations.

Transfer function measurements can be easily made from the front panel by a programmable scan of up to 11 frequencies. Setups and offsets are recalled at each frequency in the scan.

## **Analog Inputs and Outputs**

The two displays each have a user-defined output for measuring X, Y, R, R(dBm),  $\theta$ , and X-noise or Y-noise. Two user-programmable DACs provide -10.5 V to +10.5 V outputs with 1 mV resolution. These outputs may be set from the front panel or via the computer interfaces.

In addition, there are two general-purpose analog inputs. These are 16-bit ADCs which can be displayed on the front panel, read over the interface or used to ratio the input signal.

## **Internal Memory**

The SR844 has two 16,000 point memory buffers for recording (rates to 512 samples/s) the time history of each displayed measurement. Data may be transferred from the buffers using either interface. A trigger input is also provided to synchronize data recording with external events.

## **Ordering Information**

SR844 200 MHz dual phase lock-in amplifier (w/ rack mount)
SR445 Voltage preamplifier (300 MHz, 4 channel)



SR844 rear panel



#### **Signal Channel**

 $\begin{array}{lll} \mbox{Voltage input} & \mbox{Single-ended BNC} \\ \mbox{Input impedance} & \mbox{50 } \Omega \mbox{ or 1 } M\Omega + 30 \mbox{ pF} \\ \mbox{Damage threshold} & \pm 5 \mbox{ V (DC + AC)} \\ \mbox{Bandwidth} & 25 \mbox{ kHz to 200 MHz} \\ \mbox{Sensitivity} & \end{array}$ 

 $\begin{array}{lll} <1 \text{ MHz} & 100 \text{ nVrms to 1 Vrms full scale} \\ <50 \text{ MHz} & 1 \,\mu\text{Vrms to 1 Vrms full scale} \\ <200 \text{ MHz} & 10 \,\mu\text{Vrms to 1 Vrms full scale} \end{array}$ 

Gain accuracy

<50 MHz ±0.25 dB <200 MHz ±0.50 dB Gain stability 0.2 %/°C

Coherent pickup Low-noise reserve, sens. <30 mV

 $\begin{array}{lll} f<10 \ MHz & <100 \ nV \ (typ.) \\ f<50 \ MHz & <2.5 \ \mu V \ (typ.) \\ f<200 \ MHz & <25 \ \mu V \ (typ.) \end{array}$ 

Input noise (50  $\Omega$ ) 2 nV/ $\sqrt{\text{Hz}}$  (typ.), <8 nV/ $\sqrt{\text{Hz}}$  (max.) Input noise (1 M $\Omega$ ) 5 nV/ $\sqrt{\text{Hz}}$  (typ.), <8 nV/ $\sqrt{\text{Hz}}$  (max.)

Dynamic reserve up to 80 dB

### **Reference Channel**

Threshold setting Automatic, midpoint of waveform 410 s (auto-ranging, any frequency) 41 s (within same octave)

Internal reference 25 kHz to 200 MHz

Freq. resolution 3 digits

Freq. accuracy  $\pm 0.1$  in the 3<sup>rd</sup> digit Harmonic detection 2F (50 kHz to 200 MHz)

Reference outputs Phase locked to int./ext. reference 25 kHz to 200 MHz square wave

1.0 Vpp nominal into  $50 \Omega$ 25 kHz to 1.5 MHz, 0 to +5 V

Rear panel (TTL) 25 kHz to 1.5 MHz, 0 to nominal,  $\geq$ 3 V into 50  $\Omega$ 

Phase resolution 0.02°

Absolute phase error

<50 MHz <2.5° <100 MHz <5.0° <200 MHz <10.0° Rel. phase error, orthog. <2.5°

Phase noise (external) 0.005° rms at 100 MHz,

100 ms time constant

Phase drift

## **Demodulator**

Zero stability Digital displays have no zero-drift.

Analog outputs have <5 ppm/°C

drift for all dyn. reserve settings. 100 µs to 30 ks with 6, 12, 18 or

Time constants 100 µs to 30 ks with 24 dB/octave rolloff

"No Filter" mode 10 to 20 µs update rate (X and Y)

Harmonic rejection

Odd harmonics  $-9.5 \, \text{dBc} \, (20.3 \, \text{s} \, \text{ref}, -14 \, \text{dBc} \, (20.5 \, \text{s} \, \text{ref}, -14 \, \text{dBc}) = 1.5 \, \text{dBc}$ 

etc.  $(20 \log 1/n, n = 3, 5, 7...)$ 

Even harmonics < -40 dBc Sub-harmonics < -40 dBc

Spurious responses  $-10 \text{ dBc @ref } \pm 2 \times \text{IF}$ 

 $-23 \text{ dBc } @\text{ref } \pm 4 \times \text{IF}$ < -30 dBc otherwise

#### **Displays**

Channel 1

Type 4½-digit LED and 40-seg. bar graph Quantities X, R (V or dBm), X-noise, Aux In 1

Channel 2

Type  $4\frac{1}{2}$ -digit LED and 40-seg. bar graph Quantities Y,  $\theta$ , Y-noise (V or dBm), Aux In 2 Expand ×10 or ×100 for Ch1 and Ch2 Ratio X and Y ratioed with respect to

Aux In 1 or Aux In 2 before filtering and computation of R. The ratio input is normalized to 1 V and has a

dynamic range greater than 100.

Reference

Type 4½-digit LED

Quantities Ref Freq, Phase, Offsets, Aux Out,

IF Freq, Elapsed Time

## **Channel 1 and Channel 2 Outputs**

Voltage range  $\pm 10$  V full scale proportional to X,

Y or CH1, CH2 displayed quantity

Update rate

X, Y
48 to 96 kHz
R, θ, Aux inputs
X-noise, Y-noise
512 Hz

#### **Auxiliary Inputs and Outputs**

Inputs 2

Type Differential, 1 M $\Omega$ 

 $\begin{array}{ccc} \text{Range} & \pm 10 \text{ V} \\ \text{Resolution} & 0.33 \text{ mV} \\ \text{Bandwidth} & 3 \text{ kHz} \\ \text{Outputs} & 2 \\ \text{Range} & \pm 10 \text{ V} \\ \text{Resolution} & 1 \text{ mV} \end{array}$ 

Data buffers Two 16,000 point buffers. Data is

recorded at rates up to 512 Hz and is read using computer interfaces.

General

Interfaces IEEE-488.2 and RS-232 interfaces

are standard.

Power 70 W, 100/120/220/240 VAC,

50/60 Hz

Dimensions  $17" \times 5.25" \times 19.5"$  (WHD)

Weight 23 lbs.

Warranty One year parts and labor on defects

in materials and workmanship





## **About RF Lock-In Amplifiers**

The SR844 RF lock-in amplifier utilizes a combination of analog and digital techniques to obtain maximum performance over a wide frequency range. Since it is not feasible to use pure digital techniques at the SR844's maximum 200 MHz operating frequency, analog down-conversion is used to transform the signal frequency to a range suitable for DSP processing.


A block diagram of the SR844 lock-in amplifier is shown below. The RF input signal passes through adjustable RF attenuators and gain stages depending on the selected input sensitivity. The signal is then mixed with two reference signals which differ in phase by 90°.

The reference signals are generated from either the external reference input or the built-in frequency synthesizer, using a phase-lock loop (PLL) circuit. If the reference frequency were exactly at the signal frequency, the output of the mixers would be at DC. Since it is difficult to build drift and offset free amplifier chains, the SR844 chops the reference signals at a chopping frequency ( $f_c$ ) which is chosen to be fast relative to the fastest time constants, yet slow relative to the input signal frequency. The IF amplifier and filter chain can now be AC coupled, eliminating DC offset and drift problems.

Once the in-phase and quadrature IF signals have been amplified and low-pass filtered, the signals are digitized by two precision 16-bit analog-to-digital converters. The digital IF signals are ratioed (if ratio mode is selected) and digitally

low-pass filtered (with 6, 12, 18 or 24 dB/oct filter slopes) allowing the original signal amplitude and phase to be recovered.

Note that the SR844 uses a square wave mixer, not a sine wave mixer like other SRS lock-in amplifiers. This is because precision sine generation is impractical with current technology over the SR844's operating range. The effect of using a square wave mixer is that the lock-in will respond not only at the reference frequency, but at all the Fourier components of the square wave reference. Since a square wave consists of odd harmonics with amplitudes 1/3, 1/5, 1/7, etc., the SR844 will respond at odd multiples of the reference frequency as well as at the reference frequency as itself. This usually does not present a problem as long as it is understood.



SR844 block diagram



phone: (408)744-9040 www.thinkSRS.com